

Safety Data Sheet

according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 31/12/2014

Revision date: 05/13/2015

ALCOTABS

1 Identification of the Substance/mixture and of the Company/Undertaking

1.1 Product identifier

Trade name: **ALCOTABS**

Application of the substance / the preparation: Alkaline cleaner/ detergent.

1.2 Relevant identified uses of the substance or mixture and uses advised against:

No additional information available.

1.3 Details of the supplier of the Safety Data Sheet

Manufacturer/Supplier:

Alconox, Inc.
30 Glenn St., Suite 309
White Plains, NY 10603
Phone: 914-948-4040



Further information obtainable from: Product Safety Department

1.4 Emergency telephone number:

ChemTel Inc.: (800)255-3924, +1 (813)248-0585

2 Hazards Identification

2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008:



GHS07

Eye Irrit. 2A; H319: Causes serious eye irritation.

Information concerning particular hazards for human and environment:

The product has to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.

Classification system:

The classification is according to the latest editions of the EU-lists, and extended by company and literature data

2.2 Label elements

Labelling according to Regulation (EC) No 1272/2008

The product is classified and labelled according to the CLP regulation.

Hazard pictograms



GHS07

Signal word: Warning

Hazard-determining components of labelling:

Alkylbenzene Sulfonic Acid, Sodium Salt

Hazard statements:

H319: Causes serious eye irritation.

Precautionary statements:

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P264: Wash skin thoroughly after handling.

P302+352: IF ON SKIN: Wash with soap and water.

P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Safety Data Sheet

according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 31/12/2014

Revision date: 05/13/2015

ALCOTABS

Hazard description:**WHMIS classification and symbols:**

D2B - Toxic material causing other toxic effects.

**NFPA / HMIS ratings (scale 0 - 4)**

0=Minimal Hazard; 1=Slight Hazard; 2=Moderate Hazard; 3=Serious Hazard; 4=Severe Hazard; * Chronic Hazard



Health = 1
Fire = 0
Reactivity = 0

HEALTH	1	Health = 1
FIRE	0	Fire = 0
REACTIVITY	0	Reactivity = 0

2.3 Other hazards**Results of PBT and vPvB assessment**

PBT: Not applicable.

vPvB: Not applicable.

3 Composition/Information on Ingredients**3.2 Chemical characterization:** Mixture**Description:** Hazardous ingredients of mixture listed below.

Identifying Nos.	Dangerous components	Wt. %
CAS: 77-92-9 EINECS: 201-069-1	Citric acid	10 - 25%
CAS: 68081-81-2	Alkylbenzene Sulfonic Acid, Sodium Salt	2.5 - 10%
CAS: 497-19-8 EINECS: 207-838-8	Sodium carbonate	2.5 - 10%
CAS: 7722-88-5 EINECS: 231-767-1	Tetrasodium pyrophosphate	2.5 - 10%

Additional information: For the wording of the listed risk phrases refer to section 16.**4 First Aid Measures****4.1 Description of first aid measures****General information:**

Supply fresh air; consult doctor in case of complaints.

After inhalation:

Supply fresh air; consult doctor in case of complaints.

After skin contact:

Generally the product does not irritate the skin.

Immediately rinse with water.

After eye contact:

Remove contact lenses if worn.

Rinse opened eye for at least 30 minutes under running water, lifting upper and lower lids occasionally.

After swallowing:

Do not induce vomiting; call for medical help immediately.

Rinse out mouth and then drink plenty of water.

4.2 Most important symptoms and effects, both acute and delayed:

No additional information available.

4.3 Indication of any immediate medical attention and special treatment needed:

No additional information available.

Safety Data Sheet

according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 31/12/2014

Revision date: 05/13/2015

ALCOTABS

5 Firefighting Measures

5.1 Extinguishing media:

Suitable extinguishing agents:

CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

5.2 Special hazards arising from the substance or mixture:

No additional information available.

5.3 Advice for firefighters:

Protective equipment:

Wear self-contained respiratory protective device.

Wear fully protective suit.

6 Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures:

Product forms slippery surface when combined with water.

6.2 Environmental precautions: Do not allow to enter sewers/ surface or ground water.

6.3 Methods and material for containment and cleaning up:

Pick up mechanically.

Clean the affected area carefully; suitable cleaners are: Warm water

6.4 Reference to other sections:

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information

7 Handling and Storage

7.1 Precautions for safe handling:

Keep container tightly sealed.

Prevent formation of dust.

Information about fire - and explosion protection: No special measures required.

7.2 Conditions for safe storage, including any incompatibilities:

Storage:

Requirements to be met by storerooms and receptacles: No special requirements.

Information about storage in one common storage facility: Not required.

Further information about storage conditions: Protect from humidity and water.

7.3 Specific end use(s): No additional information available.

8 Exposure Controls/Personal Protection

8.1 Control parameters

Ingredients with limit values that require monitoring at the workplace:

7722-88-5 Tetrasodium pyrophosphate	
REL (USA)	5 mg/m ³
TLV (USA)	TLV withdrawn
EV (Canada)	5 mg/m ³

Additional information: The lists valid during the making were used as basis.

8.2 Exposure controls:

Personal protective equipment:

General protective and hygienic measures:

Immediately remove all soiled and contaminated clothing.

Wash hands before breaks and at the end of work.

Keep away from foodstuffs, beverages and feed.

Do not eat or drink while working.

Avoid contact with the eyes.

Safety Data Sheet

according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 31/12/2014

Revision date: 05/13/2015

ALCOTABS

Respiratory protection:

Use suitable respiratory protective device in case of insufficient ventilation.

Protection of hands:

Protective gloves

The glove material has to be impermeable and resistant to the product. Selection of the glove material should be based on the penetration time, rates of diffusion and the degradation of the glove material.

Material of gloves:

The selection of a suitable gloves does not only depend on the material, but also on the quality, and varies from manufacturer to manufacturer.

Penetration time of glove material:

The exact break through time has to be determined by the manufacturer of the protective gloves. DO NOT exceed the breakthrough time set by the Manufacturer.

For long term contact, gloves made of the following materials are considered suitable:

Butyl rubber, BR
Nitrile rubber, NBR
Natural rubber (NR)
Neoprene gloves

Eye protection:

Safety glasses
Goggles recommended during refilling

Body protection: Protective work clothing

9 Physical and Chemical Properties

9.1 Information on basic physical and chemical properties:**General Information:****Appearance:**

Form:	Solid
Color:	White
Odor:	Odourless
Odor threshold:	Not determined.
pH-value:	6.5 (- NA for Powder)

Change in condition:

Melting point/Melting range:	Undetermined.
Boiling point/Boiling range:	Undetermined.
Flash point:	Not applicable.
Flammability (solid, gaseous):	Not applicable.
Ignition temperature:	Not applicable.
Decomposition temperature:	Not determined.
Self-igniting:	Product is not selfigniting.
Danger of explosion:	Product does not present an explosion hazard.
Explosion limits:	
Lower:	Not determined.
Upper:	Not determined.
Vapor pressure at 20°C:	Not applicable.
Density:	Not applicable.
Relative density:	Not determined.
Vapor density:	Not determined.
Evaporation rate:	Not determined.
Solubility in / Miscibility with water:	Soluble.

Safety Data Sheet

according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 31/12/2014

Revision date: 05/13/2015

ALCOTABS

Segregation coefficient (n-octanol/water):	Not determined.
Viscosity:	
Dynamic:	Not applicable.
Kinematic:	Not applicable.
Solvent content:	
Organic solvents:	0.0%
Solids content:	100 %
9.2 Other information:	No additional information available

10 Stability and Reactivity

10.1 Reactivity:**10.2 Chemical stability:****Thermal decomposition / conditions to be avoided:**

No decomposition if used according to specifications.
Can decompose slowly with localized heating above 150°C.

10.3 Possibility of hazardous reactions:

Reacts with acids.
Reacts with strong alkali.
Reacts with strong oxidizing agents.

10.4 Conditions to avoid: No additional information available.**10.5 Incompatible materials:** No additional information available.**10.6 Hazardous decomposition products:**

Sulphur oxides (SO_x)
Carbon monoxide and carbon dioxide
Phosphorus compounds

11 Toxicological Information

11.1 Information on toxicological effects:**Toxicity Data:** No additional information available.**Primary irritant effect:**

On the skin: No irritant effect.
On the eye: Irritating effect.
Sensitization: No sensitizing effects known.

Additional toxicological information:

The product shows the following dangers according to the calculation method of the General EU Classification Guidelines for Preparations as issued in the latest version: Irritant.

12 Ecological Information

12.1 Toxicity:**Aquatic toxicity:** No additional information available.**12.2 Persistence and degradability:** No additional information available.**12.3 Bioaccumulative potential:** Does not accumulate in organisms.**12.4 Mobility in soil:** No additional information available.**Additional ecological information:****General notes:**

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

12.5 Results of PBT and vPvB assessment:**PBT:** Not applicable.**vPvB:** Not applicable.**12.6 Other adverse effects:** No additional information available.

Safety Data Sheet

according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 31/12/2014

Revision date: 05/13/2015

ALCOTABS

13 Disposal Considerations

13.1 Waste treatment methods:

Recommendation:

Smaller quantities can be disposed of with household waste.

Small amounts may be diluted with plenty of water and washed away. Dispose of bigger amounts in accordance with Local Authority requirements.

The surfactant used in this product complies with the biodegradability criteria as laid down in Regulation (EC) No. 648/2004 on detergents. Data to support this assertion are held at the disposal of the competent authorities of the Member States and will be made available to them, at their direct request or at the request of a detergent manufacturer.

Uncleaned packaging:

Recommendation: Disposal must be made according to official regulations.

Recommended cleansing agents: Water, together with cleansing agents, if necessary.

14 Transport Information

14.1 UN-Number:

DOT, ADR, ADN, IMDG, IATA: Not Regulated.

14.2 UN proper shipping name:

DOT, ADR, ADN, IMDG, IATA: Not Regulated.

14.3 Transport hazard class(es):

DOT, ADR, ADN, IMDG, IATA:
Class: Not Regulated.

14.4 Packing group:

DOT, ADR, ADN, IMDG, IATA: Not Regulated.

14.5 Environmental hazards:

Marine pollutant: No.

14.6 Special precautions for user: Not applicable.

14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code: Not applicable.

UN "Model Regulation": Not Regulated.

15 Regulatory Information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture:

United States (USA):

SARA:

Section 355 (extremely hazardous substances): None of the ingredient is listed.

Section 313 (Specific toxic chemical listings): None of the ingredient is listed.

TSCA (Toxic Substances Control Act): All ingredients are listed.

Proposition 65 (California):

Chemicals known to cause cancer: None of the ingredient is listed.

Chemicals known to cause reproductive toxicity for females: None of the ingredient is listed.

Chemicals known to cause reproductive toxicity for males: None of the ingredient is listed.

Chemicals known to cause developmental toxicity: None of the ingredient is listed.

Safety Data Sheet

according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 31/12/2014

Revision date: 05/13/2015

ALCOTABS

Carcinogenic Categories:

EPA (Environmental Protection Agency): None of the ingredient is listed.

TLV (Threshold Limit Value established by ACGIH): None of the ingredient is listed.

NIOSH-Ca (National Institute for Occupational Safety and Health): None of the ingredient is listed.

OSHA-Ca (Occupational Safety & Health Administration): None of the ingredient is listed.

Canadá:

Canadian Domestic Substances List (DSL): All ingredients are listed.

Canadian Ingredient Disclosure list (limit 0.1%): None of the ingredient is listed.

Canadian Ingredient Disclosure list (limit 1%):

77-92-9 citric acid

497-19-8 sodium carbonate

7722-88-5 tetrasodium pyrophosphate

15.2 Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other Information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Relevant phrases:

H319: Causes serious eye irritation.

Abbreviations and Acronyms

ADR: European Agreement concerning the International Carriage of Dangerous Goods by Road

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

ACGIH: American Conference of Governmental Industrial Hygienists

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

WHMIS: Workplace Hazardous Materials Information System (Canada)

VOC: Volatile Organic Compounds (USA, EU)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

SDS Created by:

Global Safety Management, Inc.

10006 Cross Creek Blvd

Tampa, FL, 33647

Tel: 1-844-GSM-INFO (1-844-476-4636)

Website: www.GSMSDS.com