SAFETY DATA SHEET

Ruthenium(IV) oxide hydrate

Section 1. Identification

GHS product identifier	: Ruthenium(IV) oxide hydrate
Product code	: 02596-AB, 02596-CF
Chemical name	: r úthenium oxide hydrate
Other means of identification	: Ruthenium oxide (RuO2), hydrate; Ruthenium(IV) oxide hydrate
Product type	: Solid.
Relevant identified uses of	the substance or mixture and uses advised against
Product use	: Laboratory chemicals.
Area of application	: Industrial applications, Professional applications.
Supplier's details	: SPI Supplies Division Structure Probe, Inc. 206 Garfield Ave. West Chester, PA 19380 United States Telephone: 1-(610)-436-5400 http://www.2spi.com
e-mail address of person responsible for this SDS	: SDS@2spi.com
Emergency telephone	: CHEMTREC
number (with hours of operation)	Toll Free: 1-(800)-424-9300 (USA + Canada) (24/7) International: 1-(703)-741-5970 (24/7)

Section 2. Hazards identification

OSHA/HCS status	: This material is (29 CFR 1910	s considered hazardous .1200).	by the OSHA Hazard	Communica	ation Standa	rd
Classification of the substance or mixture	: ₩315 H319 H335		0,	SINGLE EX	POSURE)	
GHS label elements						
Hazard pictograms						
Signal word	: Warning					
Hazard statements		s skin irritation. s serious eye irritation. use respiratory irritation				
Precautionary statements	,	. ,				
Date of issue/Date of revision	: 10/13/2023	Date of previous issue	: 08/31/2023	Version	:2	1/12

Section 2. Hazards identification

Prevention	 P280 - Wear protective gloves. Wear eye or face protection. P271 - Use only outdoors or in a well-ventilated area. P261 - Avoid breathing dust. P264 - Wash thoroughly after handling.
Response	 P304 - Wash thoroughly after handling. P304 + P340, P312 - IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER or doctor if you feel unwell. P362 + P364 - Take off contaminated clothing and wash it before reuse. P302 + P352 - IF ON SKIN: Wash with plenty of water. P332 + P313 - If skin irritation occurs: Get medical advice or attention. P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337 + P313 - If eye irritation persists: Get medical advice or attention.
Storage	₽405 - Store locked up. P403 + P233 - Store in a well-ventilated place. Keep container tightly closed.
Disposal	: P501 - Dispose of contents and container in accordance with all local, regional, national and international regulations.
Hazards not otherwise classified	: None known.

Section 3. Composition/information on ingredients

Substance/mixture	: Substance
Chemical name	: r uthenium oxide hydrate
Other means of identification	: Ruthenium oxide (RuO2), hydrate; Ruthenium(IV) oxide hydrate

CAS number/other identifiers

CAS number	: 32740-79-7
------------	--------------

Ingredient name	Other names	%	CAS number
ruthenium oxide hydrate	-	100	32740-79-7

Any concentration shown as a range is to protect confidentiality or is due to batch variation.

There are no additional ingredients present which, within the current knowledge of the supplier and in the concentrations applicable, are classified and hence require reporting in this section.

Occupational exposure limits, if available, are listed in Section 8.

Section 4. First aid measures

Description of necessa	ary first aid measures
Eye contact	Immediately flush eyes with plenty of water, occasionally lifting the upper and lower eyelids. Check for and remove any contact lenses. Continue to rinse for at least 10 minutes. Get medical attention.
Inhalation	: Remove victim to fresh air and keep at rest in a position comfortable for breathing. If it is suspected that fumes are still present, the rescuer should wear an appropriate mask or self-contained breathing apparatus. If not breathing, if breathing is irregular or if respiratory arrest occurs, provide artificial respiration or oxygen by trained personnel. It may be dangerous to the person providing aid to give mouth-to-mouth resuscitation. Get medical attention. If necessary, call a poison center or physician. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband.

Section 4. First aid measures

Skin contact	: Fush contaminated skin with plenty of water. Remove contaminated clothing and shoes. Continue to rinse for at least 10 minutes. Get medical attention. Wash clothing before reuse. Clean shoes thoroughly before reuse.
Ingestion	: Wash out mouth with water. Remove dentures if any. If material has been swallowed and the exposed person is conscious, give small quantities of water to drink. Stop if the exposed person feels sick as vomiting may be dangerous. Do not induce vomiting unless directed to do so by medical personnel. If vomiting occurs, the head should be kept low so that vomit does not enter the lungs. Get medical attention if adverse health effects persist or are severe. Never give anything by mouth to an unconscious person. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband.

Most important symptoms/e	ffects, acute and delayed
Potential acute health effe	<u>cts</u>
Eye contact	: 🖉 auses serious eye irritation.
Inhalation	: May cause respiratory irritation.
Skin contact	: 🖉auses skin irritation.
Ingestion	: No known significant effects or critical hazards.
<u>Over-exposure signs/symp</u>	<u>otoms</u>
Eye contact	: Adverse symptoms may include the following: pain or irritation watering redness
Inhalation	: Adverse symptoms may include the following: respiratory tract irritation coughing
Skin contact	: Adverse symptoms may include the following: irritation redness
Ingestion	: No specific data.
	lical attention and special treatment needed, if necessary
Notes to physician	 Treat symptomatically. Contact poison treatment specialist immediately if large quantities have been ingested or inhaled.
Specific treatments	: No specific treatment.
Protection of first-aiders	■ No action shall be taken involving any personal risk or without suitable training. If it is suspected that fumes are still present, the rescuer should wear an appropriate mask or self-contained breathing apparatus. It may be dangerous to the person providing aid to give mouth-to-mouth resuscitation.

See toxicological information (Section 11)

Section 5. Fire-fighting measures

Extinguishing media	
Suitable extinguishing media	: Use an extinguishing agent suitable for the surrounding fire.
Unsuitable extinguishing media	: Do not use water jet.

Date of issue/Date of revision

Section 5. Fire-fighting measures

Specific hazards arising from the chemical	: No specific fire or explosion hazard.
Hazardous thermal decomposition products	: Decomposition products may include the following materials: metal oxide/oxides
Special protective actions	: Promptly isolate the scene by removing all persons from the vicinity of the incident if there

for fire-fighters	is a fire. No action shall be taken involving any personal risk or without suitable training.
Special protective equipment for fire-fighters	: Fire-fighters should wear appropriate protective equipment and self-contained breathing apparatus (SCBA) with a full face-piece operated in positive pressure mode.

Section 6. Accidental release measures

Personal precautions, protect	tiv	e equipment and emergency procedures
For non-emergency personnel	:	No action shall be taken involving any personal risk or without suitable training. Evacuate surrounding areas. Keep unnecessary and unprotected personnel from entering. Do not touch or walk through spilled material. Provide adequate ventilation. Wear appropriate respirator when ventilation is inadequate. Put on appropriate personal protective equipment.
For emergency responders	:	If specialized clothing is required to deal with the spillage, take note of any information in Section 8 on suitable and unsuitable materials. See also the information in "For non-emergency personnel".
Environmental precautions	:	Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers. Inform the relevant authorities if the product has caused environmental pollution (sewers, waterways, soil or air).

Methods and materials for containment and cleaning up

Small spill :	Move containers from spill area. Avoid dust generation. Using a vacuum with HEPA filter will reduce dust dispersal. Place spilled material in a designated, labeled waste container. Dispose of via a licensed waste disposal contractor.
Large spill :	Move containers from spill area. Approach release from upwind. Prevent entry into sewers, water courses, basements or confined areas. Avoid dust generation. Do not dry sweep. Vacuum dust with equipment fitted with a HEPA filter and place in a closed, labeled waste container. Dispose of via a licensed waste disposal contractor. Note: see Section 1 for emergency contact information and Section 13 for waste disposal.

Section 7. Handling and storage

Precautions for safe handling

Protective measures	: Fut on appropriate personal protective equipment (see Section 8). Do not ingest. Avoid contact with eyes, skin and clothing. Use only with adequate ventilation. Wear appropriate respirator when ventilation is inadequate. Keep in the original container or an approved alternative made from a compatible material, kept tightly closed when not in use. Empty containers retain product residue and can be hazardous. Do not reuse container.
Advice on general occupational hygiene	: Eating, drinking and smoking should be prohibited in areas where this material is handled, stored and processed. Workers should wash hands and face before eating, drinking and smoking. Remove contaminated clothing and protective equipment before entering eating areas. See also Section 8 for additional information on hygiene measures.

Date of issue/Date of revision : 10/13/2023 Date of previous issue : 08/31/2023 Version : 2	Date of issue/Date of revision	: 10/13/2023 Date of previous	issue : 08/31/2023	Version : 2	4/12
---	--------------------------------	-------------------------------	--------------------	-------------	------

Section 7. Handling and storage

Conditions for safe storage,	
including any	direct sunlight in a dry, cool and well-ventilated area, away from incompatible materials
incompatibilities	(see Section 10) and food and drink. Store locked up. Keep container tightly closed and
	sealed until ready for use. Containers that have been opened must be carefully resealed
	and kept upright to prevent leakage. Do not store in unlabeled containers. Use
	appropriate containment to avoid environmental contamination. See Section 10 for
	incompatible materials before handling or use.

Section 8. Exposure controls/personal protection

Control parameters

Occupational exposure limits

Ingredient name	Exposure limits
🗖 thenium oxide hydrate	None.

Biological exposure indices

None known.

Appropriate engineering : controls	Se only with adequate ventilation. If user operations generate dust, fumes, gas, vapor or mist, use process enclosures, local exhaust ventilation or other engineering controls to keep worker exposure to airborne contaminants below any recommended or statutory limits.
Environmental exposure : controls	Emissions from ventilation or work process equipment should be checked to ensure they comply with the requirements of environmental protection legislation. In some cases, fume scrubbers, filters or engineering modifications to the process equipment will be necessary to reduce emissions to acceptable levels.
Individual protection measures	
Hygiene measures :	Wash hands, forearms and face thoroughly after handling chemical products, before eating, smoking and using the lavatory and at the end of the working period. Appropriate techniques should be used to remove potentially contaminated clothing. Wash contaminated clothing before reusing. Ensure that eyewash stations and safety showers are close to the workstation location.
Eye/face protection :	Safety eyewear complying with an approved standard should be used when a risk assessment indicates this is necessary to avoid exposure to liquid splashes, mists, gases or dusts. If contact is possible, the following protection should be worn, unless the assessment indicates a higher degree of protection: chemical splash goggles.
Skin protection	
Hand protection :	Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary. Considering the parameters specified by the glove manufacturer, check during use that the gloves are still retaining their protective properties. It should be noted that the time to breakthrough for any glove material may be different for different glove manufacturers. In the case of mixtures, consisting of several substances, the protection time of the gloves cannot be accurately estimated.
Body protection :	Personal protective equipment for the body should be selected based on the task being performed and the risks involved and should be approved by a specialist before handling this product.
Other skin protection :	Appropriate footwear and any additional skin protection measures should be selected based on the task being performed and the risks involved and should be approved by a specialist before handling this product.
Date of issue/Date of revision	: 10/13/2023 Date of previous issue : 08/31/2023 Version : 2 5/12

Section 8. Exposure controls/personal protection

Respiratory protection

: Based on the hazard and potential for exposure, select a respirator that meets the appropriate standard or certification. Respirators must be used according to a respiratory protection program to ensure proper fitting, training, and other important aspects of use.

Section 9. Physical and chemical properties

The conditions of measurement of all properties are at standard temperature and pressure unless otherwise indicated.

Appearance	
Physical state	: Solid. [Powder]
Color	: Black.
Odor	: Not available.
Odor threshold	: Not available.
рН	: Not available.
Melting point/freezing point	: Not available.
Boiling point, initial boiling point, and boiling range	: Not available.
Flash point	: Not applicable.
Flammability	: Not available.
Lower and upper explosion limit/flammability limit	: Not applicable.
Vapor pressure	: Not available.
Relative vapor density	: Not applicable.
Relative density	: 🗾
Density	: Not available.
Solubility(ies)	: Not available.
Partition coefficient: n- octanol/water	: Not available.
Auto-ignition temperature	: Not applicable.
Decomposition temperature	: Not available.
SADT	: Not available.
Viscosity	: Not applicable.
Molecular weight	: 160 g/mole
Particle characteristics	
Median particle size	: Not available.
Other information	
Physical/chemical properties comments	: No additional information

Section 10. Stability and reactivity

Reactivity	:	No specific test data related to reactivity available for this product or its ingredients.
Chemical stability	:	The product is stable.
Possibility of hazardous reactions	:	Under normal conditions of storage and use, hazardous reactions will not occur. Under normal conditions of storage and use, hazardous polymerization will not occur.

 Date of issue/Date of revision
 : 10/13/2023
 Date of previous issue

sue : 08/31/2023

Section 10. Stability and reactivity

Conditions to avoid	: Kvoid open flame and ignition sources.
Incompatible materials	: Reactive or incompatible with the following materials: organic materials: Acetone, Alcohols metals
Hazardous decomposition products	: Under normal conditions of storage and use, hazardous decomposition products should not be produced.

products

Section 11. Toxicological information

formation on toxicologica	<u>l effects</u>	
Acute toxicity		
Not available.		
Irritation/Corrosion		
Not available.		
<u>Sensitization</u>		
Not available.		
<u>Mutagenicity</u>		
Conclusion/Summary	: Not available.	
Carcinogenicity		
Conclusion/Summary	: Not available.	
Reproductive toxicity		
Conclusion/Summary	: Not available.	
Teratogenicity		
Conclusion/Summary	: Not available.	
Specific target organ toxic	<u>:ity (single exposure)</u>	

Name		Route of exposure	Target organs
ruthenium oxide hydrate	Category 3		Respiratory tract irritation

Specific target organ toxicity (repeated exposure) Not available.

Aspiration hazard

Not available.

Information on the likely	: Routes of entry anticipated: Oral, Dermal, Inhalation, Eyes.
routes of exposure	

Potential acute health effe	<u>cts</u>		
Eye contact	: 🔽auses serious eye irritation.		
Inhalation	: May cause respiratory irritation.		
Skin contact	: 🗖auses skin irritation.		
Ingestion	: No known significant effects or critical hazards.		
Date of issue/Date of revision	: 10/13/2023 Date of previous issue : 08/31/2023	Version : 2	7/12

Section 11. Toxicological information

Symptoms related to the physical, chemical and toxicological characteristics				
Eye contact		Adverse symptoms may include the following: pain or irritation watering redness		
Inhalation	:	Adverse symptoms may include the following: respiratory tract irritation coughing		
Skin contact	:	Adverse symptoms may include the following: irritation redness		
Ingestion	:	No specific data.		
	<u>ts a</u>	and also chronic effects from short and long term exposure		
Short term exposure		N1 / 11 /		
Potential immediate effects	:	Not available.		
Potential delayed effects	1	Not available.		
<u>Long term exposure</u>				
Potential immediate effects	:	Not available.		
Potential delayed effects	:	Not available.		
Potential chronic health effects				
General	1	No known significant effects or critical hazards.		
Carcinogenicity	:	No known significant effects or critical hazards.		
Mutagenicity	:	No known significant effects or critical hazards.		
Reproductive toxicity	:	No known significant effects or critical hazards.		

Numerical measures of toxicity

Acute toxicity estimates

<mark>N</mark>∕/A

Section 12. Ecological information

Toxicity

Conclusion/Summary

: Not available.

Persistence and degradability

Not available.

Bioaccumulative potential

Not available.

Mobility in soil

Date of issue/Date of revision

Section 12. Ecological information

Soil/water partition coefficient (Koc)

: Not available.

Other adverse effects

: No known significant effects or critical hazards.

Section 13. Disposal considerations

Disposal methods : The generation of waste should be avoided or minimized wherever possible. Disposal of this product, solutions and any by-products should at all times comply with the requirements of environmental protection and waste disposal legislation and any regional local authority requirements. Dispose of surplus and non-recyclable products via a licensed waste disposal contractor. Waste should not be disposed of untreated to the sewer unless fully compliant with the requirements of all authorities with jurisdiction. Waste packaging should be recycled. Incineration or landfill should only be considered when recycling is not feasible. This material and its container must be disposed of in a safe way. Care should be taken when handling emptied containers that have not been cleaned or rinsed out. Empty containers or liners may retain some product residues. Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers.

Section 14. Transport information

	DOT Classification	IMDG	ΙΑΤΑ		
UN number	Not regulated.	Not regulated.	Not regulated.		
UN proper shipping name					
Transport hazard class(es)					
Packing group					
Environmental hazards	No.	No.	No.		

Additional information

Special precautions for user : Transport within user's premises: always transport in closed containers that are upright and secure. Ensure that persons transporting the product know what to do in the event of an accident or spillage.

Transport in bulk according : Not available. to IMO instruments

Date of issue/Date of revision

Section 15. Regulatory information

0				
U.S. Federal regulations		ι,	npt/Partial exemption: Not determined	
	Un	ited States invente	ory (TSCA 8b): This material is active or exempted.	
Clean Air Act Section 112 (b) Hazardous Air Pollutants (HAPs)	: No	Not listed		
Clean Air Act Section 602 Class I Substances	: No	Not listed		
Clean Air Act Section 602 Class II Substances	: Na	t listed		
DEA List I Chemicals (Precursor Chemicals)	: No	Not listed		
DEA List II Chemicals (Essential Chemicals)	: No	Not listed		
SARA 302/304				
Composition/information	on ingr	edients		
No products were found.				
SARA 304 RQ	: No	t applicable.		
SARA 311/312				
Classification : SKIN IRRITATION - Category 2 EYE IRRITATION - Category 2A SPECIFIC TARGET ORGAN TOXICITY (SINGLE EXPOSURE) (Respiratory tract irritation) - Category 3				
Composition/information on ingredients				
Name		%	Classification	
		100		

Name	%	Classification
rdthenium oxide hydrate		SKIN IRRITATION - Category 2 EYE IRRITATION - Category 2A SPECIFIC TARGET ORGAN TOXICITY (SINGLE EXPOSURE) (Respiratory tract irritation) - Category 3

SARA 313

Not applicable.

State regulations

Massachusetts	: This material is not listed.
New York	: This material is not listed.
New Jersey	: This material is listed.

Pennsylvania : This material is listed.

California Prop. 65

This product does not require a Safe Harbor warning under California Prop. 65.

International regulations

Chemical Weapon Convention List Schedules I. II & III Chemicals Not listed.

Montreal Protocol

Not listed.

Stockholm Convention on Persistent Organic Pollutants

Not listed.

```
Date of issue/Date of revision
```

: 10/13/2023 Date of previous issue

: 08/31/2023

Section 15. Regulatory information

Rotterdam Convention on Prior Informed Consent (PIC) Not listed.

UNECE Aarhus Protocol on POPs and Heavy Metals Not listed.

Section 16. Other information

Hazardous Material Information System (U.S.A.)



Caution: HMIS® ratings are based on a 0-4 rating scale, with 0 representing minimal hazards or risks, and 4 representing significant hazards or risks. Although HMIS® ratings and the associated label are not required on SDSs or products leaving a facility under 29 CFR 1910.1200, the preparer may choose to provide them. HMIS® ratings are to be used with a fully implemented HMIS® program. HMIS® is a registered trademark and service mark of the American Coatings Association, Inc.

The customer is responsible for determining the PPE code for this material. For more information on HMIS® Personal Protective Equipment (PPE) codes, consult the HMIS® Implementation Manual.

National Fire Protection Association (U.S.A.)

Health Plammability Health Plammability/Reactivity Special hazards

Procedure used to derive the classification

	Classification	Justification
KIN IRRITATION - Catego EYE IRRITATION - Categor SPECIFIC TARGET ORGAN irritation) - Category 3		Expert judgment Expert judgment Expert judgment
History		· ·
Date of issue/Date of revision	: 10/13/2023	
Date of previous issue	: 08/31/2023	
Version	: 2	
Prepared by	: Sphera Solutions	
Key to abbreviations	: ATE = Acute Toxicity Estimate AMP = Acceptable maximum peak above the accep 8-hr shift BCF = Bioconcentration Factor GHS = Globally Harmonized System of Classificatio IATA = International Air Transport Association IBC = International Air Transport Association IBC = International Maritime Dangerous Goods LogPow = logarithm of the octanol/water partition co MARPOL = International Convention for the Prevent modified by the Protocol of 1978. ("Marpol" = marine N/A = Not available UN = United Nations	n and Labelling of Chemicals efficient ion of Pollution From Ships, 1973 as
Date of issue/Date of revision	: 10/13/2023 Date of previous issue : 08/31/2023	Version : 2 11/12

Section 16. Other information

References

: HCS (U.S.A.) - Hazard Communication Standard International transport regulations

Indicates information that has changed from previously issued version.

Notice to reader

To the best of our knowledge, the information contained herein is accurate. However, neither the above-named supplier, nor any of its subsidiaries, assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.