


TECHNICAL NOTE

TAPERED CRUCIBLES (Type T)

SPI Supplies
206 Garfield Avenue,
West Chester, PA 19380, USA

PRODUCT DIMENSION INFORMATION

SPI-Glas™ 22 only occurs inside on the crucible, not on the exterior of the side walls for the crucible.



Available Products:

Volume ml	D mm	B mm	H mm	Wall Thickness, mm	SPI #
0.07	6	4	11	1	4238GCCT-AB
17	31	26	43	2.5	4239GCCT-AB
18	24	15	87	2	4240GCCT-AB
37	51	42	30	3	4241GCCT-AB
64	45	44	58	3	4242GCCT-AB
105	59	49	59	3	4243GCCT-AB
190	45	41	168	3	4244GCCT-AB
227	73	68	85	3	4245GCCT-AB
320	73	70	100	3	4246GCCT-AB
370	69	52	133	2	4247GCCT-AB
550	85	75	130	2	4248GCCT-AB
3060	132	122	273	2	4249GCCT-AB

Revised: EER 1/16