

# Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Keylajet®

## I Identification of the substance/mixture and of the supplier

### I.1 GHS Product identifier

Trade Name: Keylajet®

Product number: 2405, 2415, 2455

### I.2 Application of the substance / the mixture: Cleaning material/Detergent

### I.3 Details of the supplier of the Safety Data Sheet

**Manufacturer:**

Alconox Inc.  
30 Glenn St  
White Plains, NY 10603  
(914) 948-4040

**Supplier:****Emergency telephone number:**

ChemTel Inc  
North America: 1-888-255-3924  
International: +1 813-248-0573

## 2 Hazards identification

### 2.1 Classification of the substance or mixture:

In compliance with EC regulation No. 1272, 29CFR1910/1200 and GHS and requirements.

**Hazard-determining components of labeling:**

Sodium hydroxide  
Potassium hydroxide  
Tetrasodium EDTA  
Tetrapotassium pyrophosphate

### 2.2 Label elements:

Skin corrosion, category 1A.  
Serious eye damage, category 1.

**Hazard pictograms:****Signal word:** Danger**Hazard statements:**

H314 Causes severe skin burns and eye damage.  
H318 Causes serious eye damage.

**Precautionary statements:**

P264 Wash skin thoroughly after handling.  
P280 Wear protective gloves/protective clothing/eye protection/face protection.  
P301+P330+P312 IF SWALLOWED: Rinse mouth. Call a POISON CENTER or doctor/physician if you feel unwell.  
P302+P352+P312 If on skin: Wash with soap and water. Call a poison center or doctor/physician if you feel unwell.  
P304+P341 If inhaled: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.  
P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing.  
P337+P313 If eye irritation persists get medical advice/attention.  
P501 Dispose of contents and container as instructed in Section 13.

## Safety Data Sheet

**Effective date:** 11 May 2020

**Revision:** 11 May 2020

**Trade Name:** Keylajet®

**Additional information:** None.

### Hazard description

**Hazards Not Otherwise Classified (HNOC):** May cause surfaces to become slippery. Use caution in areas of foot traffic if on floors.

### Information concerning particular hazards for humans and environment:

The product has to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.

### Classification system:

The classification is according to EC regulation No. 1272, 29CFR1910/1200 and GHS and requirements, and extended by company and literature data. The classification is in accordance with the latest editions of international substances lists, and is supplemented by information from technical literature and by information provided by the company.

## 3 Composition/information on ingredients

**3.1 Chemical characterization:** None

**3.2 Description:** None

**3.3 Hazardous components (percentages by weight)**

Identification	Chemical Name	Classification	Wt. %
<b>CAS number:</b> 1310-73-2	Sodium Hydroxide	Acute Tox. 4; H303 Skin Corr. 1A; H314	1-10
<b>CAS number:</b> 1310-58-3	Potassium Hydroxide	Acute Tox. 4; H303 Skin Corr. 1A; H314	1-20
<b>CAS number:</b> 7320-34-5	Tetrapotassium Pyrophosphate	Eye irritation, Cat. 2B	3-10
<b>CAS number:</b> 64-02-8	Tetrasodium EDTA	Acute Tox. 4; H303 Eye Dam. 1; H318	1-5

**3.4 Additional Information:** None.

## 4 First aid measures

**4.1 Description of first aid measures**

**General information:** None.

### After inhalation:

Move exposed individual to fresh air.  
Loosen clothing as necessary and position individual in a comfortable position.  
Maintain an unobstructed airway.  
Immediately call a POISON CONTROL CENTER or seek medical attention.

### After skin contact:

Immediately remove all contaminated clothing.  
Wash affected area with soap and water.  
Immediately call a POISON CONTROL CENTER or seek medical attention.

### After eye contact:

Rinse/flush exposed eye(s) gently using water for 15-20 minutes.  
Remove contact lens(es) if able to do so during rinsing.  
Immediately call a POISON CONTROL CENTER or seek medical attention.

## Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Keylajet®

**After swallowing:**

Immediately call a POISON CONTROL CENTER or seek medical attention.  
Do not induce vomiting.  
Rinse mouth and then drink plenty of water.

**4.2 Most important symptoms and effects, both acute and delayed**

None

**4.3 Indication of any immediate medical attention and special treatment needed:**

No additional information.

### 5 Firefighting measures

**5.1 Extinguishing media**

**Suitable extinguishing agents:**

Use appropriate fire suppression agents for adjacent combustible materials or sources of ignition.

**For safety reasons unsuitable extinguishing agents:** None

**5.2 Special hazards arising from the substance or mixture :**

Thermal decomposition can lead to release of irritating gases and vapors.

**5.3 Advice for firefighters**

**Protective equipment:**

Wear protective eye wear, gloves and clothing.  
Refer to Section 8.

**5.4 Additional information:**

Avoid inhaling gases, fumes, dust, mist, vapor and aerosols.  
Avoid contact with skin, eyes and clothing.

### 6 Accidental release measures

**6.1 Personal precautions, protective equipment and emergency procedures:**

Ensure adequate ventilation.  
Ensure air handling systems are operational.

**6.2 Environmental precautions:**

Should not be released into the environment.  
Prevent from reaching drains, sewer or waterway.

**6.3 Methods and material for containment and cleaning up:**

Wear protective eye wear, gloves and clothing.

**6.4 Reference to other sections:** None

### 7 Handling and storage

**7.1 Precautions for safe handling:**

Avoid vapor, if atomized.  
Do not eat, drink, smoke or use personal products when handling chemical substances.

**7.2 Conditions for safe storage, including any incompatibilities:**

Store in a cool, well-ventilated area.

**7.3 Specific end use(s):**

No additional information.

## Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Keylajet®

### 8 Exposure controls/personal protection



#### 8.1 Control parameters:

1310-73-2, Sodium hydroxide, NIOSH C 2.0 mg/m<sup>3</sup>.1310-58-3, Sodium hydroxide, ACGIH C 2.0 mg/m<sup>3</sup>. This material will not become airborne under normal operating conditions.

#### 8.2 Exposure controls

##### Appropriate engineering controls:

Emergency eye wash fountains and safety showers should be available in the immediate vicinity of use or handling.

##### Respiratory protection:

Not needed under normal conditions.

##### Protection of skin:

Select glove material impermeable and resistant to the substance – nitrile or equivalent recommended.

##### Eye protection:

Safety goggles or glasses, or appropriate eye protection.

##### General hygienic measures:

Wash hands before breaks and at the end of work.

Avoid contact with skin, eyes and clothing.

### 9 Physical and chemical properties

<b>Appearance (physical state, color):</b>	Clear pale yellow liquid	<b>Explosion limit lower:</b> <b>Explosion limit upper:</b>	Not determined or not available. Not determined or not available.
<b>Odor:</b>	Not determined or not available.	<b>Vapor pressure at 20°C:</b>	Not determined or not available.
<b>Odor threshold:</b>	Not determined or not available.	<b>Vapor density:</b>	Not determined or not available.
<b>pH-value:</b>	13+ (1% aqueous solution)	<b>Relative density:</b>	Not determined or not available.
<b>Melting/Freezing point:</b>	Not determined or not available.	<b>Solubilities:</b>	Not determined or not available.
<b>Boiling point/Boiling range:</b>	Not determined or not available.	<b>Partition coefficient (n-octanol/water):</b>	Not determined or not available.
<b>Flash point (closed cup):</b>	Not determined or not available.	<b>Auto/Self-ignition temperature:</b>	Not determined or not available.
<b>Evaporation rate:</b>	Not determined or not available.	<b>Decomposition temperature:</b>	Not determined or not available.
<b>Flammability (solid, gaseous):</b>	Not determined or not available.	<b>Viscosity:</b>	a. Kinematic: Not determined or not available. b. Dynamic: Not determined or not available.
<b>Density at 20°C:</b>	Not determined or not available.		

## Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Keylajet®

### 10 Stability and reactivity

- 10.1 Reactivity:** Not determined or not available.
- 10.2 Chemical stability:** Not determined or not available.
- 10.3 Possibility hazardous reactions:** Not determined or not available.
- 10.4 Conditions to avoid:** Not determined or not available.
- 10.5 Incompatible materials:** Avoid mixing with strong organic and mineral acids, reducing agents.
- 10.6 Hazardous decomposition products:** Not determined or not available.

### 11 Toxicological information

#### 11.1 Information on toxicological effects:

**Acute Toxicity:****Oral:**

: LD50 &gt;500 mg/kg oral rat - product.

**Chronic Toxicity:** No additional information.**Skin corrosion/irritation:**

Sodium hydroxide: Severe skin irritation - 24h.

**Serious eye damage/irritation:**

Sodium hydroxide: Corrosive to eyes.

**Respiratory or skin sensitization:** No additional information.**Carcinogenicity:** No additional information.**IARC (International Agency for Research on Cancer):** None of the ingredients are listed.**NTP (National Toxicology Program):** None of the ingredients are listed.**Germ cell mutagenicity:** No additional information.**Reproductive toxicity:** No additional information.**STOT-single and repeated exposure:** No additional information.**Additional toxicological information:** No additional information.

### 12 Ecological information

#### 12.1 Toxicity:

Sodium hydroxide: Fish, LC50 - Gambusia affinis (Mosquito fish) - 125 mg/l - 96 h.

Potassium hydroxide: Fish, LC50 - Gambusia affinis (Mosquito fish) - 80 mg/l - 96 h.

**12.2 Persistence and degradability:** No additional information.**12.3 Bioaccumulative potential:** No additional information.**12.4 Mobility in soil:** No additional information.**General notes:** No additional information.

#### 12.5 Results of PBT and vPvB assessment:

**PBT:** No additional information.**vPvB:** No additional information.**12.6 Other adverse effects:** No additional information.

### 13 Disposal considerations

#### 13.1 Waste treatment methods (consult local, regional and national authorities for proper disposal)

**Relevant Information:**

It is the responsibility of the waste generator to properly characterize all waste materials according to applicable regulatory entities. (US 40CFR262.11).

## Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Keylajet®

**I4 Transport information**

<b>I4.1 UN Number:</b> ADR, ADN, DOT, IMDG, IATA	1760						
<b>I4.2 UN Proper shipping name:</b> ADR, ADN, DOT, IMDG, IATA	Corrosive Liquid, N.O.S. (Potassium Hydroxide and Sodium Hydroxide Solution).						
<b>I4.3 Transport hazard classes:</b> ADR, ADN, DOT, IMDG, IATA	<table> <tr> <td><b>Class:</b></td> <td>8</td> </tr> <tr> <td><b>Label:</b></td> <td>8 Corrosive substances</td> </tr> <tr> <td><b>LTD.QTY:</b></td> <td>None</td> </tr> </table>	<b>Class:</b>	8	<b>Label:</b>	8 Corrosive substances	<b>LTD.QTY:</b>	None
<b>Class:</b>	8						
<b>Label:</b>	8 Corrosive substances						
<b>LTD.QTY:</b>	None						
<p><b>US DOT</b></p> <p><b>Limited Quantity Exception:</b> None</p> <table> <tr> <td> <p><b>Bulk:</b></p> <p><b>RQ (if applicable):</b> None</p> <p><b>Proper shipping Name:</b> Corrosive Liquid, N.O.S. (Potassium Hydroxide and Sodium Hydroxide Solution).</p> <p><b>Hazard Class:</b> 8</p> <p><b>Packing Group:</b> II.</p> <p><b>Marine Pollutant (if applicable):</b> No</p> <p><b>Comments:</b> None</p> </td> <td> <p><b>Non Bulk:</b></p> <p><b>RQ (if applicable):</b> None</p> <p><b>Proper shipping Name:</b> Corrosive Liquid, N.O.S. (Potassium Hydroxide and Sodium Hydroxide Solution).</p> <p><b>Hazard Class:</b> 8</p> <p><b>Packing Group:</b> II.</p> <p><b>Marine Pollutant (if applicable):</b> No</p> <p><b>Comments:</b> None</p> </td> </tr> </table>		<p><b>Bulk:</b></p> <p><b>RQ (if applicable):</b> None</p> <p><b>Proper shipping Name:</b> Corrosive Liquid, N.O.S. (Potassium Hydroxide and Sodium Hydroxide Solution).</p> <p><b>Hazard Class:</b> 8</p> <p><b>Packing Group:</b> II.</p> <p><b>Marine Pollutant (if applicable):</b> No</p> <p><b>Comments:</b> None</p>	<p><b>Non Bulk:</b></p> <p><b>RQ (if applicable):</b> None</p> <p><b>Proper shipping Name:</b> Corrosive Liquid, N.O.S. (Potassium Hydroxide and Sodium Hydroxide Solution).</p> <p><b>Hazard Class:</b> 8</p> <p><b>Packing Group:</b> II.</p> <p><b>Marine Pollutant (if applicable):</b> No</p> <p><b>Comments:</b> None</p>				
<p><b>Bulk:</b></p> <p><b>RQ (if applicable):</b> None</p> <p><b>Proper shipping Name:</b> Corrosive Liquid, N.O.S. (Potassium Hydroxide and Sodium Hydroxide Solution).</p> <p><b>Hazard Class:</b> 8</p> <p><b>Packing Group:</b> II.</p> <p><b>Marine Pollutant (if applicable):</b> No</p> <p><b>Comments:</b> None</p>	<p><b>Non Bulk:</b></p> <p><b>RQ (if applicable):</b> None</p> <p><b>Proper shipping Name:</b> Corrosive Liquid, N.O.S. (Potassium Hydroxide and Sodium Hydroxide Solution).</p> <p><b>Hazard Class:</b> 8</p> <p><b>Packing Group:</b> II.</p> <p><b>Marine Pollutant (if applicable):</b> No</p> <p><b>Comments:</b> None</p>						
<b>I4.4 Packing group:</b> ADR, ADN, DOT, IMDG, IATA	II.						
<b>I4.5 Environmental hazards:</b>	None						
<b>I4.6 Special precautions for user:</b>	None						
<b>Danger code (Kemler):</b>	None						
<b>EMS number:</b>	None						
<b>Segregation groups:</b>	None						
<b>I4.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code:</b> Not applicable.							
<b>I4.8 Transport/Additional information:</b>							
<b>Transport category:</b>	None						
<b>Tunnel restriction code:</b>	None						
<b>UN "Model Regulation":</b>	UN1760, Corrosive Liquid, N.O.S. (Potassium Hydroxide and Sodium Hydroxide Solution), 8, PG II						

## Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Keylajet®

### 15 Regulatory information

#### 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture.

##### North American

###### SARA

**Section 313 (specific toxic chemical listings):** None of the ingredients are listed.

**Section 302 (extremely hazardous substances):** None of the ingredients are listed.

###### CERCLA (Comprehensive Environmental Response, Clean up and Liability Act) Reportable

**Spill Quantity:** None of the ingredients are listed.

###### TSCA (Toxic Substances Control Act):

**Inventory:** All ingredients are listed as active.

**Rules and Orders:** Not applicable.

###### Proposition 65 (California):

**Chemicals known to cause cancer:** None of the ingredients are listed.

**Chemicals known to cause reproductive toxicity for females:** None of the ingredients are listed.

**Chemicals known to cause reproductive toxicity for males:** None of the ingredients are listed.

**Chemicals known to cause developmental toxicity:** None of the ingredients are listed.

###### Canadian

###### Canadian Domestic Substances List (DSL):

All ingredients are listed.

##### EU

**REACH Article 57 (SVHC):** None of the ingredients are listed.

**Germany MAK:** Not classified.

##### Asia Pacific

###### Australia

**Australian Inventory of Chemical Substances (AICS):** All ingredients are listed.

###### China

**Inventory of Existing Chemical Substances in China (IECSC):** All ingredients are listed.

###### Japan

**Inventory of Existing and New Chemical Substances (ENCS):** All ingredients are listed.

###### Korea

**Existing Chemicals List (ECL):** All ingredients are listed.

###### New Zealand

**New Zealand Inventory of Chemicals (NZOIC):** All ingredients are listed.

###### Taiwan

**Taiwan Chemical Substance Inventory (TSCI):** All ingredients are listed.

## Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Keylajet®

### 16 Other information

**Abbreviations and Acronyms:** None

#### Summary of Phrases

**Hazard statements:**

H314 Causes severe skin burns and eye damage.

H318 Causes serious eye damage.

**Precautionary statements:**

P264 Wash skin thoroughly after handling.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P330+P312 IF SWALLOWED: Rinse mouth. Call a POISON CENTER or doctor/physician if you feel unwell.

P302+P352+P312 If on skin: Wash with soap and water. Call a poison center or doctor/physician if you feel unwell.

P304+P341 If inhaled: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.

P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing.

P337+P313 If eye irritation persists get medical advice/attention.

P501 Dispose of contents and container as instructed in Section 13.

#### Manufacturer Statement:

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

**NFPA:** 2-0-0

**HMIS:** 2-0-0