

Safety Data Sheet

according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 11/20/2014

Revision: 05/12/2015

ALCOJET

1 Identification of the Substance/mixture and of the Company/Undertaking

1.1 Product identifier

Trade name: **ALCOJET**

Application of the substance / the preparation: Cleaning material/ Detergent

1.2 Relevant identified uses of the substance or mixture and uses advised against:

No additional information available.

1.3 Details of the supplier of the Safety Data Sheet

Manufacturer/Supplier:

Alconox, Inc.
30 Glenn St., Suite 309
White Plains, NY 10603
Phone: 914-948-4040



Further information obtainable from: Product Safety Department

1.4 Emergency telephone number:

ChemTel Inc.: (800)255-3924, +1 (813)248-0585

2 Hazards Identification

2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008:



GHS07

*Eye Irrit. 2A; H319: Causes eye irritation.**Skin Irrit. 2; H315: Causes skin irritation.**STOT SE 3; H335: May cause respiratory irritation.*

Information concerning particular hazards for human and environment:

The product has to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.

Classification system:

The classification is according to the latest editions of the EU-lists, and extended by company and literature data

2.2 Label elements

Labelling according to Regulation (EC) No 1272/2008:

The product is classified and labelled according to the CLP regulation.

Hazard pictograms:



GHS07

Signal word: Warning

Hazard-determining components of labelling:

Disodium Metasilicate

Hazard statements:

H315: Causes skin irritation.

H319: Causes serious eye irritation.

H335: May cause respiratory irritation.

Safety Data Sheet

according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 11/20/2014

Revision: 05/12/2015


ALCOJET

Precautionary statements:

- P261: Avoid breathing dust/fume/gas/mist/vapors/spray.
 P280: Wear protective gloves/protective clothing/eye protection/face protection.
 P273: Avoid release to the environment.
 P264: Wash thoroughly after handling.
 P271: Use only outdoors or in a well-ventilated area.
 P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
 P362: Take off contaminated clothing and wash before reuse.
 P304+P340: IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
 P332+P313: If skin irritation occurs: Get medical advice/attention.
 P302+P352: IF ON SKIN: Wash with plenty of soap and water.
 P403+P233: Store in a well-ventilated place. Keep container tightly closed.
 P501: Dispose of contents/container in accordance with local/regional/national/international regulations.

Other Hazard description:**WHMIS-classification and symbols:**

D2B - Toxic material causing other toxic effects

**NFPA ratings (scale 0 - 4)**

Health = 1
 Fire = 0
 Reactivity = 0

HMIS-ratings (scale 0 - 4)

| | | |
|------------|---|--|
| HEALTH | 1 | |
| FIRE | 0 | |
| REACTIVITY | 0 | |

Health = 1
 Fire = 0
 Reactivity = 0

2.3 Other hazards**Results of PBT and vPvB assessment**

PBT: Not applicable.

vPvB: Not applicable.

3 Composition/Information on Ingredients**3.2 Chemical characterization:** Mixture**Description:** Hazardous ingredients of mixture listed below.

| Identifying Nos. | Description | Wt. % |
|--------------------------------------|------------------------------|-----------|
| CAS: 1344-09-8 EINECS: 215-687-4 | Silicic Acid, Sodium Salt | 10 - 25% |
| CAS: 6834-92-0 EINECS: 229-912-9 | Disodium Metasilicate | 2.5 - 10% |
| CAS: 7722-88-5 EINECS: 231-767-1 | Tetrasodium Pyrophosphate | 2.5 -10% |
| CAS: 51580-86-0 EINECS: 220-767-7 | Troclosene Sodium, Dihydrate | ≤ 2.5% |
| CAS: N/A | Proprietary (non-classified) | 25 - 40% |

Additional information: For the wording of the listed risk phrases refer to section 16.

Safety Data Sheet
according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 11/20/2014

Revision: 05/12/2015

ALCOJET**4 First Aid Measures****4.1 Description of first aid measures****General information:**

Symptoms of poisoning may even occur after several hours; therefore medical observation for at least 48 hours after the accident.

After inhalation:

Supply fresh air; consult doctor in case of complaints.

In case of unconsciousness place patient stably in side position for transportation.

After skin contact:

Immediately wash with water and soap and rinse thoroughly. If skin irritation continues, consult a doctor.

After eye contact:

Remove contact lenses if worn.

Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

After swallowing:

Do not induce vomiting; call for medical help immediately.

4.2 Most important symptoms and effects, both acute and delayed:

No additional information available.

4.3 Indication of any immediate medical attention and special treatment needed:

No additional information available.

5 Firefighting Measures**5.1 Extinguishing media:****Suitable extinguishing agents:**

CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

5.2 Special hazards arising from the substance or mixture:

No additional information available.

5.3 Advice for firefighters:**Protective equipment:**

Wear self-contained respiratory protective device.

Wear fully protective suit.

6 Accidental Release Measures**6.1 Personal precautions, protective equipment and emergency procedures:**

Product forms slippery surface when combined with water.

Avoid formation of dust.

6.2 Environmental precautions:

Do not allow product to reach sewage system or any water course.

6.3 Methods and material for containment and cleaning up:

Pick up mechanically.

Clean the affected area carefully; suitable cleaners are: Warm water

Dispose contaminated material as waste according to item 13.

Ensure adequate ventilation.

6.4 Reference to other sections:

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information

Safety Data Sheet
 according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
 GHS

Effective date: 11/20/2014

Revision: 05/12/2015

ALCOJET

7 Handling and Storage**7.1 Precautions for safe handling:**

- Ensure good ventilation/exhaustion at the workplace.
- Keep receptacles tightly sealed.
- Prevent formation of dust.

Information about fire - and explosion protection: No special measures required.

7.2 Conditions for safe storage, including any incompatibilities:**Storage:**

Requirements to be met by storerooms and receptacles: No special requirements.

Information about storage in one common storage facility: Store away from foodstuffs.

Further information about storage conditions:

- Keep container tightly sealed.
- Protect from humidity and water.
- Store in dry conditions.

7.3 Specific end use(s): No additional information available.

8 Exposure Controls/Personal Protection**8.1 Control parameters**

Ingredients with limit values that require monitoring at the workplace: Not required.

| | |
|-----------|---------------------------|
| 7722-88-5 | tetrasodium pyrophosphate |
|-----------|---------------------------|

| | |
|-------------|---------------------|
| REL (USA) | 5 mg/m ³ |
| TLV (USA) | TLV withdrawn |
| EV (Canada) | 5 mg/m ³ |

Additional information: The lists valid during the making were used as basis.

8.2 Exposure controls:**Personal protective equipment:****General protective and hygienic measures:**

- Keep away from foodstuffs, beverages and feed.
- Immediately remove all soiled and contaminated clothing.
- Wash hands before breaks and at the end of work.
- Avoid contact with the eyes and skin.

Respiratory protection:

- In case of brief exposure or low pollution use respiratory filter device.
- In case of intensive or longer exposure use self-contained respiratory protective device.

Protection of hands:

Protective gloves

The glove material has to be impermeable and resistant to the product. Selection of the glove material should be based on the penetration time, rates of diffusion and the degradation of the glove material.

Material of gloves:

The selection of a suitable gloves does not only depend on the material, but also on the quality, and varies from manufacturer to manufacturer.

Penetration time of glove material:

The exact break through time has to be determined by the manufacturer of the protective gloves. DO NOT exceed the breakthrough time set by the Manufacturer.

Safety Data Sheet

according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 11/20/2014

Revision: 05/12/2015

ALCOJET

For long term contact, gloves made of the following materials are considered suitable:

Butyl rubber, BR
Nitrile rubber, NBR
Natural rubber (NR)
Neoprene gloves

Eye protection:



Safety glasses

Body protection: Protective work clothing

9 Physical and Chemical Properties

9.1 Information on basic physical and chemical properties:

General Information:

Appearance:

| | |
|------------------------|-----------------|
| Form: | Powder |
| Color: | White |
| Odor: | Characteristic |
| Odor threshold: | Not determined. |
| pH-value: | NA for Powder. |

Change in condition:

| | |
|-------------------------------------|-----------------|
| Melting point/Melting range: | Not determined. |
| Boiling point/Boiling range: | Not determined. |

Flash point: Not applicable.

Flammability (solid, gaseous): Not applicable.

Ignition temperature: Not determined.

Decomposition temperature: Not determined.

Self-igniting: Product is not selfigniting.

Danger of explosion: Product does not present an explosion hazard.

Explosion limits:

| | |
|---------------|-----------------|
| Lower: | Not determined. |
| Upper: | Not determined. |

Vapor pressure: Not determined.

Density: Not determined.

Relative density: Not determined.

Vapor density: Not applicable.

Evaporation rate: Not applicable.

Solubility in / Miscibility with water: Fully miscible.

Segregation coefficient (n-octanol/water): Not determined.

Viscosity:

| | |
|-------------------|-----------------|
| Dynamic: | Not determined. |
| Kinematic: | Not determined. |

Solvent content:

Organic solvents: 0.0 %

Solids content: 100 %

9.2 Other information: No additional information available.

Safety Data Sheet

according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 11/20/2014

Revision: 05/12/2015

ALCOJET

10 Stability and Reactivity

10.1 Reactivity:**10.2 Chemical stability:****Thermal decomposition / conditions to be avoided:**

No decomposition if used according to specifications.

Can decompose slowly with localized heating above 150°C.

10.3 Possibility of hazardous reactions:

Contact with acids releases irritant gases.

Reacts with acids.

10.4 Conditions to avoid: No additional information available.**10.5 Incompatible materials:** No additional information available.**10.6 Hazardous decomposition products:**

Carbon monoxide and carbon dioxide.

Sulphur dioxide.

Nitrogen oxides.

Sulphur trioxide (SO₃) or SO₃-mist.

Phosphorus compounds.

11 Toxicological Information

11.1 Information on toxicological effects:**Acute toxicity:**

| LD/LC50 values relevant for classification: | | |
|---|------------------|------------------|
| 6834-92-0 Disodium Metasilicate | | |
| Oral | LD ₅₀ | 1280 mg/kg (rat) |
| 51580-86-0 Troclosen Sodium, Dihydrate | | |
| Oral | LD ₅₀ | 1400 mg/kg (rat) |

Primary irritant effect:**On the skin:** Irritating to skin and mucous membranes.**On the eye:** Strong irritant with the danger of severe eye injury.**Sensitization:** No sensitizing effects known.**Additional toxicological information:**

The product shows the following dangers according to the calculation method of the General EU Classification Guidelines for Preparations as issued in the latest version:

Harmful

Irritant

12 Ecological Information

12.1 Toxicity:**Aquatic toxicity:** No additional information available.**12.2 Persistence and degradability:** No additional information available.**12.3 Bioaccumulative potential:** Not worth-mentioning accumulating in organisms.**12.4 Mobility in soil:** No additional information available.**Ecotoxicological effects:****Remark:** Harmful to fish**Additional ecological information:****General notes:**

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

Rinse off of bigger amounts into drains or the aquatic environment may lead to increased pH-values. A high pH-value harms aquatic organisms. In the dilution of the use-level the pH-value is considerably reduced, so that after the use of the product the aqueous waste, emptied into drains, is only low water-dangerous.

Safety Data Sheet

according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 11/20/2014

Revision: 05/12/2015

ALCOJET

12.5 Results of PBT and vPvB assessment:

PBT: Not applicable.

vPvB: Not applicable.

12.6 Other adverse effects: No additional information available.

13 Disposal Considerations

13.1 Waste treatment methods:**Recommendation:**

Smaller quantities can be disposed of with household waste.

Can be disposed of with household garbage after consulting with the waste disposal facility operator and the pertinent authorities and adhering to the necessary technical regulations.

Small amounts may be diluted with plenty of water and washed away. Dispose of bigger amounts in accordance with Local Authority requirements.

The surfactant used in this product complies with the biodegradability criteria as laid down in Regulation (EC) No. 648/2004 on detergents. Data to support this assertion are held at the disposal of the competent authorities of the Member States and will be made available to them, at their direct request or at the request of a detergent manufacturer.

Uncleaned packaging:**Recommendation:** Disposal must be made according to official regulations.**Recommended cleansing agents:** Water, together with cleansing agents, if necessary.

14 Transport Information

14.1 UN-Number:

DOT, ADR, ADN, IMDG, IATA: Not Regulated

14.2 UN proper shipping name:

DOT, ADR, IMDG, IATA: Not Regulated

14.3 Transport hazard class(es):

DOT, ADR, IMDG, IATA:

Class: Not Regulated

Label: -

14.4 Packing group:

DOT, ADR, IMDG, IATA: Not Regulated

14.5 Environmental hazards:

Marine pollutant: No

14.6 Special precautions for user:

Not applicable.

14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code: Not applicable.

UN "Model Regulation": Not Regulated

15 Regulatory Information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture:

United States (USA):

SARA:

Section 355 (extremely hazardous substances): None of the ingredient is listed.**Section 313 (Specific toxic chemical listings):** None of the ingredient is listed.**TSCA (Toxic Substances Control Act):** All ingredients are listed.

Safety Data Sheet

according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 11/20/2014

Revision: 05/12/2015

ALCOJET

Proposition 65 (California):

Chemicals known to cause cancer: None of the ingredient is listed.

Chemicals known to cause reproductive toxicity for females: None of the ingredient is listed.

Chemicals known to cause reproductive toxicity for males: None of the ingredient is listed.

Chemicals known to cause developmental toxicity: None of the ingredient is listed.

Carcinogenic Categories:

EPA (Environmental Protection Agency): None of the ingredient is listed.

TLV (Threshold Limit Value established by ACGIH): None of the ingredient is listed.

NIOSH-Ca (National Institute for Occupational Safety and Health): None of the ingredient is listed.

OSHA-Ca (Occupational Safety & Health Administration): None of the ingredient is listed.

Canadá:

Canadian Domestic Substances List (DSL): All ingredients are listed.

Canadian Ingredient Disclosure list (limit 0.1%): None of the ingredient is listed.

Canadian Ingredient Disclosure list (limit 1%):

6834-92-0 Disodium Metasilicate

7722-88-5 Tetrasodium Pyrophosphate

15.2 Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other Information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Relevant phrases:

H315: Causes skin irritation.

H319: Causes serious eye irritation.

H335: May cause respiratory irritation.

Abbreviations and Acronyms:

ADR: European Agreement concerning the International Carriage of Dangerous Goods by Road.

IMDG: International Maritime Code for Dangerous Goods.

DOT: US Department of Transportation.

IATA: International Air Transport Association.

GHS: Globally Harmonized System of Classification and Labelling of Chemicals.

ACGIH: American Conference of Governmental Industrial Hygienists.

NFPA: National Fire Protection Association (USA).

HMIS: Hazardous Materials Identification System (USA).

WHMIS: Work-place Hazardous Materials Information System (Canada).

VOC: Volatile Organic Compounds (USA, EU).

LC50: Lethal concentration, 50 percent.

LD50: Lethal dose, 50 percent.

Safety Data Sheet
according to 1907/2006/EC (REACH), 1272/2008/EC (CLP), and
GHS

Effective date: 11/20/2014

Revision: 05/12/2015

ALCOJET

SDS Created by:

Global Safety Management, Inc.
10006 Cross Creek Blvd
Tampa, FL, 33647
Tel: 1-844-GSM-INFO (1-844-476-4636)
Website: www.GSMSDS.com