

Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Detergent 8®

I Identification of the substance/mixture and of the supplier

I.1 GHS Product identifier

Trade Name: Detergent 8®

Product number: 1701, 1701-I, 1705, 1715, 1755

I.2 Application of the substance / the mixture: Cleaning material/Detergent

I.3 Details of the supplier of the Safety Data Sheet

Manufacturer	Supplier
Alconox Inc. 30 Glenn St White Plains, NY 10603 (914) 948-4040	

Emergency telephone number:

ChemTel Inc

North America: 1-888-255-3924

International: +1 813-248-0573

2 Hazards identification

2.1 Classification of the substance or mixture:

In compliance with EC regulation No. 1272, 29CFR1910/1200 and GHS requirements.

Hazard-determining components of labeling:

2-butoxyethanol

1-aminopropan-2-ol

2.2 Label elements:

Skin corrosion, category 1C.

Serious eye damage, category 1.

Hazard pictograms:



Signal word: Danger

Hazard statements:

H314 Causes severe skin burns and eye damage.

Precautionary statements:

P264 Wash skin thoroughly after handling.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P330+P312 IF SWALLOWED: Rinse mouth. Call a POISON CENTER or doctor/physician if you feel unwell.

P302+P352+P312 If on skin: Wash with soap and water. Call a poison center or doctor/physician if you feel unwell.

P304+P341 If inhaled: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.

Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Detergent 8®

P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing.

P337+P313 If eye irritation persists get medical advice/attention.

P501 Dispose of contents and container as instructed in Section 13.

Additional information: None.

Hazard description

Hazards Not Otherwise Classified (HNOC): None

Information concerning particular hazards for humans and environment:

The product has to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.

Classification system:

The classification is according to EC regulation No. 1272/2008, 29CFR1910/1200 and GHS Rev. 3 and amendments, and extended by company and literature data. The classification is in accordance with the latest editions of international substances lists and is supplemented by information from technical literature and by information provided by the company.

3 Composition/information on ingredients

3.1 Chemical characterization: None

3.2 Description: None

3.3 Hazardous components (percentages by weight)

Identification	Chemical Name	Classification	Wt. %
CAS number: 78-96-6	1-aminopropan-2-ol	Acute Tox. 4; H313 Skin Corr. 1B; H314 Eye Dam. 1; H318	25-50
CAS number: 111-76-2	2-butoxyethanol	Acute Tox. 4; H303 Acute Tox. 4; H313 Acute Tox. 4; H333 Skin Irrit. 2; H315 Eye Irrit. 2; H319	2.5-10

3.4 Additional Information: None.

4 First aid measures

4.1 Description of first aid measures

General information: None.

After inhalation:

Move exposed individual to fresh air.

Loosen clothing as necessary and position individual in a comfortable position.

Maintain an unobstructed airway.

Immediately call a POISON CONTROL CENTER or seek medical attention.

After skin contact:

Immediately remove all contaminated clothing.

Wash affected area with soap and water.

Immediately call a POISON CONTROL CENTER or seek medical attention.

After eye contact:

Rinse/flush exposed eye(s) gently using water for 15-20 minutes.

Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Detergent 8®

Remove contact lens(es) if able to do so during rinsing.
Immediately call a POISON CONTROL CENTER or seek medical attention.

After swallowing:

Immediately call a POISON CONTROL CENTER or seek medical attention.
Do not induce vomiting.
Rinse mouth and then drink plenty of water.

4.2 Most important symptoms and effects, both acute and delayed

None

4.3 Indication of any immediate medical attention and special treatment needed:

No additional information.

5 Firefighting measures

5.1 Extinguishing media

Suitable extinguishing agents:

Use appropriate fire suppression agents for adjacent combustible materials or sources of ignition.

For safety reasons unsuitable extinguishing agents: None

5.2 Special hazards arising from the substance or mixture:

Thermal decomposition can lead to release of ammonia containing irritating gases and vapors.

5.3 Advice for firefighters

Protective equipment:

Wear protective eye wear, gloves and clothing.
Refer to Section 8.

5.4 Additional information:

Avoid inhaling gases, fumes, dust, mist, vapor and aerosols.
Avoid contact with skin, eyes and clothing.

6 Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures :

Ensure adequate ventilation.
Ensure air handling systems are operational.

6.2 Environmental precautions:

Should not be released into the environment.
Prevent from reaching drains, sewer or waterway.

6.3 Methods and material for containment and cleaning up :

Wear protective eye wear, gloves and clothing.

6.4 Reference to other sections: None

7 Handling and storage

7.1 Precautions for safe handling:

Avoid breathing mist or vapor.
Do not eat, drink, smoke or use personal products when handling chemical substances.

7.2 Conditions for safe storage, including any incompatibilities:

Store in a cool, well-ventilated area.

7.3 Specific end use(s):

No additional information.

Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Detergent 8®

8 Exposure controls/personal protection



8.1 Control parameters:

- 111-76-2, 2-butoxyethanol, ACGIH TWA 20.0 ppm.
- 111-76-2, 2-butoxyethanol, OSHA TWA 50.0 ppm 240.0 mg/m³.
- 111-76-2, 2-butoxyethanol, NIOSH TWA 5.0 ppm 24.0 mg/m³.

8.2 Exposure controls

Appropriate engineering controls:

Emergency eye wash fountains and safety showers should be available in the immediate vicinity of use or handling.

Respiratory protection:

Not needed under normal conditions.

Protection of skin:

Select glove material impermeable and resistant to the substance.

Eye protection:

Safety goggles or glasses, or appropriate eye protection.

General hygienic measures:

Wash hands before breaks and at the end of work.

Avoid contact with skin, eyes and clothing.

9 Physical and chemical properties

Appearance (physical state, color):	Clear liquid	Explosion limit lower: Explosion limit upper:	Not determined or not available. Not determined or not available.
Odor:	Not determined or not available.	Vapor pressure at 20°C:	Not determined or not available.
Odor threshold:	Not determined or not available.	Vapor density:	Not determined or not available.
pH-value:	11 (1% aqueous solution)	Relative density:	Not determined or not available.
Melting/Freezing point:	Not determined or not available.	Solubilities:	Not determined or not available.
Boiling point/Boiling range:	Not determined or not available.	Partition coefficient (n-octanol/water):	Not determined or not available.
Flash point (closed cup):	Not determined or not available.	Auto/Self-ignition temperature:	Not determined or not available.
Evaporation rate:	Not determined or not available.	Decomposition temperature:	Not determined or not available.
Flammability (solid, gaseous):	Not determined or not available.	Viscosity:	a. Kinematic: Not determined or not available. b. Dynamic: Not determined or not available.
Density at 20°C:	Not determined or not available.		

10 Stability and reactivity

10.1 Reactivity: Not determined or not available.

Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Detergent 8®

10.2 Chemical stability: Not determined or not available.

10.3 Possibility hazardous reactions: Not determined or not available.

10.4 Conditions to avoid: Not determined or not available.

10.5 Incompatible materials: Not determined or not available.

10.6 Hazardous decomposition products: Not determined or not available.

11 Toxicological information

11.1 Information on toxicological effects:

Acute Toxicity:

Oral:

: LD50 > 500 mg/kg oral rat. - Product.

Chronic Toxicity: No additional information.

Skin corrosion/irritation:

2-butoxyethanol: skin irritation - 20h.

1-aminopropan-2-ol: mild skin irritation - 24h.

Serious eye damage/irritation:

2-butoxyethanol: eye irritation - 24h.

1-aminopropan-2-ol: Severe eye irritation - 24h.

Respiratory or skin sensitization:

2-butoxyethanol: Skin irritation - 20 h.

2-butoxyethanol: Eye irritation - 24 h.

Carcinogenicity: No additional information.

IARC (International Agency for Research on Cancer):

Group 3 - Not classifiable as to its carcinogenicity to humans: 2-butoxyethanol.

NTP (National Toxicology Program): None of the ingredients are listed.

Germ cell mutagenicity: No additional information.

Reproductive toxicity:

2-butoxyethanol: Overexposure may cause reproductive disorder(s) based on tests with laboratory animals.

STOT-single and repeated exposure: No additional information.

Additional toxicological information: No additional information.

12 Ecological information

12.1 Toxicity:

1-aminopropan-2-ol: Fish, LC50 - Carassius auratus (goldfish) - 210 mg/l - 96 h.

Dipropylene Glycol Monomethyl Ether: Aquatic invertebrates, EC50 - Daphnia magna (Water flea) - 1,919 mg/l - 48 h.

2-butoxyethanol: Aquatic Plants, Pseudokirchneriella subcapitata (green algae) - 911mg/l - 72 h. 2-

butoxyethanol: Aquatic invertebrates, Daphnia magna (Water flea) - 1,550 mg/l - 48 h.

2-butoxyethanol: Fish, Oncorhynchus mykiss (rainbow trout) - 1,474 mg/l - 96 h.

12.2 Persistence and degradability: No additional information.

12.3 Bioaccumulative potential: No additional information.

Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Detergent 8®
12.4 Mobility in soil: No additional information.

General notes: No additional information.

12.5 Results of PBT and vPvB assessment:
PBT: No additional information.

vPvB: No additional information.



12.6 Other adverse effects: No additional information.

13 Disposal considerations

13.1 Waste treatment methods (consult local, regional and national authorities for proper disposal)
Relevant Information:

It is the responsibility of the waste generator to properly characterize all waste materials according to applicable regulatory entities. (US 40CFR262.11).

14 Transport information

14.1 UN Number:	1760
ADR, ADN, DOT, IMDG, IATA	
14.2 UN Proper shipping name:	Corrosive Liquid, N.O.S. (Monoisopropanolamine Solution)
ADR, ADN, DOT, IMDG, IATA	
14.3 Transport hazard classes:	
ADR, ADN, DOT, IMDG, IATA	
Class:	8
Label:	8 Corrosive substances
LTD.QTY:	None
US DOT	
Limited Quantity Exception:	None
Bulk:	Non Bulk:
RQ (if applicable): None	RQ (if applicable): None
Proper shipping Name: Corrosive Liquid, N.O.S. (Monoisopropanolamine Solution).	Proper shipping Name: Corrosive Liquid, N.O.S. (Monoisopropanolamine Solution).
Hazard Class: 8	Hazard Class: 8
Packing Group: III.	Packing Group: III.
Marine Pollutant (if applicable): No	Marine Pollutant (if applicable): No
Comments: None	Comments: None
	
14.4 Packing group:	III.
ADR, ADN, DOT, IMDG, IATA	
14.5 Environmental hazards:	None
14.6 Special precautions for user:	None
Danger code (Kemler):	None

Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Detergent 8®

EMS number: None

Segregation groups: None

14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code: Not applicable.

14.8 Transport/Additional information:

Transport category: None

Tunnel restriction code: None

UN "Model Regulation": UN1760, Corrosive Liquid, N.O.S.
(Monoisopropanolamine Solution), 8, P G III

15 Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture.

North American

SARA

Section 313 (specific toxic chemical listings):

111-76-2 2-butoxyethanol.

Section 302 (extremely hazardous substances): None of the ingredients are listed.

CERCLA (Comprehensive Environmental Response, Clean up and Liability Act) Reportable

Spill Quantity: None of the ingredients are listed.

TSCA (Toxic Substances Control Act):

Inventory: All ingredients are listed as active.

Rules and Orders: Not applicable.

Proposition 65 (California):

Chemicals known to cause cancer: None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for females: None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for males: None of the ingredients are listed.

Chemicals known to cause developmental toxicity: None of the ingredients are listed.

Canadian

Canadian Domestic Substances List (DSL):

All ingredients are listed.

EU

REACH Article 57 (SVHC): None of the ingredients are listed.

Germany MAK: Not classified.

Asia Pacific

Safety Data Sheet

Effective date: 11 May 2020

Revision: 11 May 2020

Trade Name: Detergent 8®

Australia

Australian Inventory of Chemical Substances (AICS): All ingredients are listed.

China

Inventory of Existing Chemical Substances in China (IECSC): All ingredients are listed.

Japan

Inventory of Existing and New Chemical Substances (ENCS): All ingredients are listed.

Korea

Existing Chemicals List (ECL): All ingredients are listed.

New Zealand

New Zealand Inventory of Chemicals (NZOIC): All ingredients are listed.

Philippines

Philippine Inventory of Chemicals and Chemical Substances (PICCS): All ingredients are listed.

Taiwan

Taiwan Chemical Substance Inventory (TSCI): All ingredients are listed.

16 Other information

Abbreviations and Acronyms: None

Summary of Phrases

Hazard statements:

H314 Causes severe skin burns and eye damage.

Precautionary statements:

P264 Wash skin thoroughly after handling.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P330+P312 IF SWALLOWED: Rinse mouth. Call a POISON CENTER or doctor/physician if you feel unwell.

P302+P352+P312 If on skin: Wash with soap and water. Call a poison center or doctor/physician if you feel unwell.

P304+P341 If inhaled: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.

P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing.

P337+P313 If eye irritation persists get medical advice/attention.

P501 Dispose of contents and container as instructed in Section 13.

Manufacturer Statement:

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

NFPA: 2-0-0

HMIS: 2-0-0