

Safety Data Sheet

Complies with EC no. 1907/2006

Issue Date:08/15/2018

Section 1: Chemical Product and Company Identification

Cat#: 15923C

Page 1 of 3

Part Name: INITIATOR C (USE CAT.#HXL012) LOWICRYL K4M

Supplier: Polysciences, Inc.
400 Valley Road
Warrington, PA 18976 USA
MSDS Telephone #215-343-6484 Emergency only #215-378-4526

Identified uses: Laboratory use, manufacture of substances

Section 2: Hazards Identification

Hazard Overview

Harmful if inhaled or swallowed.

GHS Classification

Acute Toxicity Oral Cat 4 , Acute Toxicity Inhalation Cat 4

Signal word: **Warning**



Hazard and Precautionary Statements

H302 + H332	Harmful if swallowed or if inhaled
P280A	Wear protective gloves
P285	In case of inadequate ventilation wear respiratory protection.
P301D	IF SWALLOWED, induce vomiting as directed by medical personnel.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P305B	IF IN EYES, Separate eyelids with finger tips.
P313	Get medical advice/attention.
P335	Brush off loose particles from skin.
P351	Rinse cautiously with water for several minutes.
P501	Dispose of contents/container to proper waste area in accordance with institutional practices and local, state or federal regulations.

NFPA Rating

Hazard Ratings:

These ratings are Polysciences' Inc. own assesments of the properties of the material using the ANSI/NFPA 704 Standard. Additional information can be found by consulting in the NFPA published ratings lists (List 325 and List 49).

If no data is listed the information is not available.

Health	Flammability	Reactivity
2	1	0

Section 3: Composition/ Information on Ingredients

Note: Items listed with a CASRN... number have no CAS# available.

Item#	Name	EINECS	CAS#	% in product
1	Acetophenone, 2-methoxy-2-phenyl	222-538-7	0003524627	95-100

Section 4: First Aid Measures

Brush off the skin.

Contact medical personnel.

Flush eyes with flowing water for at least 15 minutes.

If swallowed, induce vomiting as directed by medical personnel.

Separate eyelids with finger tips.

Wash skin with deluge of water for at least 15 minutes.

Section 5: Fire Fighting Measures

Flash point, deg F.: no data

Method: no data

UEL: no data LEL: no data

Autoignition temperature, deg. F.: no data

Flammability Classification: nap

Flame Propagation Rate: no data

Hazardous Combustion Products: none

Section 6: Accidental Release Measures

Any information listed below is to be considered in addition to internal guidelines for isolation of spill, containment of spill, removal of ignition sources from immediate area, and collection for disposal of spill by trained, properly protected clean up personnel.

No special measures are indicated.

Section 7: Handling and Storage

Store at room temp

Section 8: Exposure Controls/ Personal Protection

OSHA (ACGIH) Exposure Limits

		TWA		STEL		CEILING	
		ppm	mg/	ppm	mg/m3	ppm	mg/m3
CAS#: 0003524627	IDLH: NE						
OSHA		NE	NE	NE	NE	NE	NE
ACGIH		NE	NE	NE	NE	NE	NE

The use of eye protection in the form of safety glasses with side shields and the use of skin protection for hands in the form of gloves are considered minimum and non-discretionary in work places and laboratories. Any recommended personal protection equipment or environmental equipment is to be considered as additional to safety glasses and gloves.

Use latex or equivalent gloves.

Use process enclosures, local exhaust ventilation, or other engineering controls.

Chemical-resistant gloves should be worn whenever this material is handled. The glove material has to be impermeable and resistant to the product. Gloves should be removed and replaced immediately if there is any indication of degradation or chemical breakthrough.

Rinse and remove gloves immediately after use. Wash hands with soap and water. All glove recommendations presume that the risk of exposure is through splash and not intentional immersion of the hands into the product.

Since glove permeation data does not exist for this material, no recommendation for the glove material can be given for the product.

Permeation data must be obtained from the glove manufacturer to determine if the glove is suitable for the task.

Section 9: Physical and Chemical Properties

Formula:	C15H14O2	vapor pressure:	no data
Formula Weight:	226.3	vapor density:	no data
boiling point:	no data	Specific gravity:	0
melting point:	122	ph:	no data
solubility:	no data	appearance:	crystalline powder

Section 10: Stability and Reactivity

Chemical Stabilit stable

Conditions to Avoid: heat

Incompatibility with other materials: oxidizers

Hazardous Decomposition Products: none

Hazardous Polymerization: will not occur

Section 11: Toxicological Information

Acute Data: LD50 oral rat 7560 mg/kg

Subchronic data: no data

Section 12: Ecological Information

no data

Section 13: Disposal Considerations

The following chart lists the status of the chemical and its components in reference to 40 CFR Part 261.33. If the product is listed by code number the substance may be subject to special federal and state disposal regulations. If no codes are listed the material must be disposed in compliance with all Federal, State and Local Regulations.

CAS#	Waste Code	Regulated Name
0003524627	not listed	not listed

Section 14: Transportation Data

Proper Shipping Name NOT REGULATED

Chemical Name

UN

Class

PG

Section 15: Regulatory Information

All components of this product are on the TSCA public inventory.

All components of this product are on the TSCA public inventory.

Prop 65 - Column A identifies those items which are known to the State of California to cause cancer. Column B identified items which are known to the State of California to cause reproductive toxicity.

CAS#	Column A	Column B
0003524627	no	no

State Regulatory Information :If a CAS# is listed below this material is subject to the listed state right-to-know requirements.

CAS#
0003524627 not listed

SARA Toxic Release Chemicals(as defined in Section 313 of SARA Title III)

This list identifies the toxic chemicals, including their de minimis concentrations for which reporting is required under Section 313 of the Emergency Planning and Community Right-to-Know Act (EPCRA). The list is also referred to as the Toxics Release Inventory (TRI) List.

CAS#	Regulated name	de minimis conc. %	Rep. Thres.
0003524627	not listed	not listed	not listed

SARA Extremely Hazardous Substances and TPQs

This list includes hazardous chemicals as defined in 29 CFR 1910.1200(c); and extremely hazardous substances regulated under Section 302 of SARA Title III with their TPQs (in pounds), as listed in 40 CFR 355, Appendices A and B.

CAS#	Regulated name	TPQ (pounds)	EHS-RQ(pounds)
0003524627	not listed	not listed	not listed

CERCLA
The hazardous substances, and their reportable quantities (RQs) are listed in the federal regulations at 40 CFR Part 302, Table 302.4. Release of a CERCLA hazardous substance in an amount equal to or greater than its RQ, in any 24-hour period, must be reported to the National Response Center at (800) 424-8802.

CAS#	Regulated name	RQ (pounds)
0003524627	Not listed	Not listed

Section 16: Other Information

POLYSCIENCES, INC. provides the information contained herein in good faith but makes no representation as to its comprehensiveness or accuracy. Individuals receiving this information must exercise their independent judgment in determining its appropriateness for a particular purpose. POLYSCIENCES, INC. makes no representations or warranties, either expressed or implied of merchantability, fitness for particular purposes with respect to the information set forth herein or to which the information refers. Accordingly, POLYSCIENCES, INC. will not be responsible for damages resulting from the use of or reliance upon this information.

END OF MSDS