

CARGILLE LIQUIDS SHELF LIFE



SPI Supplies
206 Garfield Avenue,
West Chester, PA 19380, USA

Cargille™ Liquids – Shelf Life

Cargille immersion oils and liquids have a suggested shelf life. Liquids should be stored away from sunlight when not in use. The shelf life of an opened bottle should be no more than half that of an unopened bottle.

Cargille Product Shelf Life

Product Name	Shelf Life Unopened
Immersion Oils Type A, B, NVH, OVH, HF, FF, 37, 300	10 Years
Immersion Oils Types: LDF, 37LDF	5 Years
Refractive Index Liquids Series AAA, AA, A, B	5 Years
R.I. Series M and Gem Refractometer	18 Months
R.I. Liquids Series E: Indices 1.500-1.570	2 1/2 Years
Series E: Indices 1.575-1.640	5 Years
Meltpoints Code: 53, 5870, 25761	20 Years
Optical Gels Codes: 0607, 0608	10 Years

EER

Date: 12-20