


TECHNICAL NOTE

TAPERED CRUCIBLES (Type K)

SPI Supplies
206 Garfield Avenue,
West Chester, PA 19380, USA

PRODUCT DIMENSION INFORMATION

The SPI-Glas™ Glassy Carbon Crucibles, Tapered Type "K", are cylindrical in shape and with a flat bottom. The SPI-Glas™ Brand Glassy Carbon Crucibles should never be used for the melting of chromium-cobalt-molybdenum alloys. Similarly, they should not be used for melting anything that easily carbides, for example silicon or tungsten.



Model	Volume ml	D mm	B mm	H mm	Wall Thickness, mm	SPI #
KAG 8	8	31	17	27	3	4250GCCK-AB
KAG 17	17	35	18	38	3	4251GCCK-AB
KAG 35	35	44	23	48	3	4252GCCK-AB
KAG 67	67	52	26	56	3	4253GCCK-AB
KAG 116	116	70	34	52	3	4254GCCK-AB
KAG 143	143	64	33	73	3	4255GCCK-AB
KAG 1220	1220	130	62	128	3	4256GCCK-AB
KAG 3000	3000	130	86	176	3	4257GCCK-AB

*D = Outside Diameter at top

*B = Outside Diameter at bottom

*H = Height

Revised: EER 1/16